

Atlas Heavy Duty Lift Chair – CH4100

The Atlas is a wide seat, high weight capacity infinite position lift chair with high quality Chenille upholstery. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease movement of the chair.

With a X-Wide/Average size rating, this chair will accommodate people of wider widths and average heights.



- Three position recline
- Four locking swivel casters included



The spring assisted footrest on all models will flip up if it comes in contact with an obstruction when lowering to minimize injury or damage.



Supplied locking swivel casters can be installed on all models, except Earth and Mars which have 3" casters. Ideal for users who may want to move their chair occasionally or rental units.

FEATURES:

- Three position function
- Chenille upholstery
- Four casters (2 rear casters braked)
- Lounger style that fully supports the legs (no gap between footrest and seat)
- Removable back for easy handling
- Easy to use hand control
- Superb stability in full lift position
- Excellent side to side stability
- Engineered bushings and washers at pivot points
- provide quiet smooth operation and durability
- Direct drive activation synchronizes right and left sides for easy opening
- Low voltage motor, 24 volt, UL approved
- Battery backup
- The Atlas has two Dewert motors that work together to run this chair like a single motor lift chair.
- Maximum weight capacity is 210kgs (462 lbs).



MODEL NAMES	ATLAS
Model numbers:	CH4100
Weight Capacity:	462 lbs (210 kg)
Seat Width:	25.2" (64 cm)
Seat Depth:	21.25" (54 cm)
Back Height:	27.6" (70 cm)
Seat Height:	19.3" (49 cm)
Overall Width:	42.1" (108 cm)
Overall Depth:	3.3" (82 cm)
Overall Height (down):	44" (115 cm)
Overall Height (raised):	58.5" (150 cm)
Max. Vertical Lift:	n/a
Box 1	108 x 82 x 66 cm
Box 2	95 x 83 x 33 cm
Net weight	60 kg (box 1) 12 kg (box 2)
Gross Weight	65 kg (box 1) 15 kg (box 2)