

# JUPITER LIFT CHAIR



## Jupiter Two Motor CH4002 (Blue) CH4003 (Brown) CH4019 (Dark Brown) CH4024 (Grey)

The Jupiter lift chair is an infinite position lift chair available with brown, dark brown or blue upholstery. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease repositioning the chair.

With a Medium/Average size rating, this chair will accommodate people of medium width and average height.



- Three position recline
- Four locking swivel casters included



The spring assisted footrest on all models will flip up if it comes in contact with an obstruction when lowering to minimize injury or damage.



Supplied locking swivel casters can be installed on all models, except Earth and Mars which have 3" casters. Ideal for users who may want to move their chair occasionally or rental units.

### FEATURES:

- Lounger style that fully supports the legs (no gap between footrest & seat)
- Removable back for easy handling
- Easy to use hand control
- UL / CUL Approved
- Maximum seat elevation with superb stability in full lift position
- Wide construction provides greater side to side stability
- Top steel base for durability & cross supports
- Engineered bushings & washers at pivot points provide quiet smooth operation & durability
- Direct drive activation synchronizes right & left sides for easy opening
- Low voltage motor, 24 volt
- Battery Back Up



MODEL NAMES	JUPITER
Model numbers:	CH4002 (Blue) CH4003 (Brown) CH4019 (Dark Brown) CH4024 (Grey)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	19.5" (49.5 cm)
Seat Depth:	20.3" (51.6 cm)
Back Height:	28.1" (71.4 cm)
Seat Height:	19.5" (49.5 cm)
Overall Width:	35.1" (89.2 cm)
Overall Depth:	73.3" (186.2 cm)
Overall Height (down):	43.7" (111 cm)
Overall Height (raised):	58.9" (149.6 cm)
Max. Vertical Lift:	n/a
Box 1	85cm x 94cm x 69cm
Box 2	79cm x 83cm x 35cm
Net weight	50 kg (box1) 8.5 kg (box 2)
Gross Weight	54 kg (box1) 10.5 kg (box 2)